INVOICE

DATE	02/08/2023
INVOICE #	502023
PO No #	871/RSMU/UMKEU/VII/2023

BILL TO:

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RAKA RADITYA KUSUMA INDRA

Maintenance Check List


Operating Microscope Type : **M501**
Customer : **RS UNDAAN SURABAYA**
Location : **KAWAR OPERASI**

Serial No. Stand : **151104004**
Maintenance Date : **29/07/2023**

Step	Description	Required Result	O.K.	Fails
1	General			
1.1	lock stand with two locking knobs.	Optics carrier movement is locked for up/down and forward/backward direction.	✓	
1.2	Loosen footbrake to "free" position	The stand can be moved easily.	✓	
1.3	Lock foot break to	The stand stays stable, the footbrake lever is locked.	✓	
2	Illumination			
2.1	Remove main lamp 1	Warning, the lamp housing, the lamp is hot! Check button appear, press check button, the error check main lamp 1 is displayed. Build in the main lamp 1 and commit the error messages.	✓	
2.2	Switch to main lamp 1 and remove main lamp 2	Warning, the lamp housing, the lamp is hot! Check button appear, press check button, the error check main lamp 2 is displayed. Build in the main lamp 2 and commit the error messages.	✓	
2.3	Observe field of illumination; set diaphragm of optics carrier to maximum	Homogeneous illumination in the complete field.	✓	
2.4	Press the buttons for Light intensity control to the maximum for main lamp 1 and afterwards for main lamp 2	Intensity increases to the maximum of 100% for main lamp 1 and afterwards the same procedure for main lamp 2. Please note the Ottoflex lamp must be set to 0% intensity.	✓	
2.5	Press the button for Light intensity control to the maximum for Ottoflex	Intensity increases to the maximum of 100% for Ottoflex. Please note the Main Lamp 1 or main lamp 2 must be set to 0% intensity.	NA	
2.6	Set the main lamp (1 or 2) to 50% and the Ottoflex lamp as well to 50%.	It is possible to set both lamp to the same value.	NA	
2.7	Observe field of illumination; set diaphragm of optics carrier to maximum	Homogeneous illumination in the complete field.	✓	
2.8	Activate the depth enhancer (button on the control units touch panel)	Depth enhancer in activated (hear a click and see in the inlet of optics carrier the field of view is smaller)	NA	
3	Brakes (only for electromagnetic brakes)			
3.1	Press side button of right handgrip.	Brakes are released for rotation movements of the microscope (break of pivot 5 standard setting)	NA	
3.2	Press side button of left handgrip.	Brakes are released for rotation movements of the microscope (break of pivot 5 standard setting)	NA	
3.3	Press rear button of left handgrip.	All brakes are released; microscope can be moved in all axes	NA	
4	Optic carrier			
4.1	Alignment/ adjustment		✓	
4.2	Contamination of the path of observation		✓	
4.3	Dove tail interface: Fitting Contamination Loose screws		✓	
4.4	Continuity of the zoom's movement (min and max speed)		✓	
4.5	Does the zoom move throughout the whole range		✓	

User :

Technician :

RS MATA UNDAAN
29/7/23

SEM


Maintenance Check List

 Operating Microscope Type : **M822**
 Customer : **RS UNDAAN SURABAYA**
 Location : **LABOR OPERASI**

 Serial No. Stand : **021 215001**
 Maintenance Date : **29/07/2023**



Step	Description	Required Result	O.K.	Fails
1	General			
1.1	Lock stand with two locking knobs.	Optics carrier movement is locked for up/down and forward/backward direction	✓	
1.2	Loosen footbrake to "free" position	The stand can be moved easily.	✓	
1.3	Lock foot brake to	The stand stays stable, the footbrake lever is locked.	✓	
2	Illumination			
2.1	Remove main lamp 1	Warning, the lamp housing, the lamp is hot! Check button appear, press check button, the error check main lamp 1 is displayed. Build in the main lamp 1 and commit the error messages.	✓	
2.2	Switch to main lamp 1 and remove main lamp 2	Warning, the lamp housing, the lamp is hot! Check button appear, press check button, the error check main lamp 2 is displayed. Build in the main lamp 2 and commit the error messages.	✓	
2.3	Observe field of illumination; set diaphragm of optics carrier to maximum	Homogeneous illumination in the complete field.	✓	
2.4	Press the buttons for Light intensity control to the maximum for main lamp 1 and afterwards for main lamp 2	Intensity increases to the maximum of 100% for main lamp 1 and afterwards the same procedure for main lamp 2. Please note the Ottoflex lamp must be set to 0% intensity.	✓	
2.5	Press the button for Light intensity control to the maximum for Ottoflex	Intensity increases to the maximum of 100% for Ottoflex. Please note the Main Lamp 1 or main lamp 2 must be set to 0% intensity.	✓	
2.6	Set the main lamp (1 or 2) to 50% and the Ottoflex lamp as well to 50%.	It is possible to set both lamp to the same value.	✓	
2.7	Observe field of illumination; set diaphragm of optics carrier to maximum	Homogeneous illumination in the complete field.	✓	
2.8	Activate the depth enhancer (button on the control units touch panel)	Depth enhancer is activated (hear a click and see in the inlet of optics carrier the field of view is smaller)	✓	
3	Brakes (only for electromagnetic brakes)			
3.1	Press side button of right handgrip.	Brakes are released for rotation movements of the microscope (break of pivot 5 standard setting)	NA	
3.2	Press side button of left handgrip.	Brakes are released for rotation movements of the microscope (break of pivot 5 standard setting)	NA	
3.3	Press rear button of left handgrip.	All brakes are released; microscope can be moved in all axes	NA	
4	Optic carrier			
4.1	Alignment/ adjustment		✓	
4.2	Contamination of the path of observation		✓	
4.3	Dove tail interface: Fitting Contamination Loose screws		✓	
4.4	Continuity of the zoom's movement (min and max speed)		✓	
4.5	Does the zoom move throughout the whole range		✓	

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